

OLYMPUS®



EVIS EXERA COLONOVIDEOSCOPE
OLYMPUS PCF TYPE 160AL/I

SLIMMER TUBE WITH VARIABLE STIFFNESS




Breakthrough innoflex™ design lets you adjust insertion tube flexibility to match the internal contours of the colon for more precise control while narrow 11.5 mm diameter insertion tube assures smooth insertion

innoflex™



Olympus business areas



OLYMPUS®
FOCUS ON LIFE

<http://www.olympus.com>

OLYMPUS OPTICAL CO., LTD.
San-Ei Building, 22-2, Nishi Shinjuku 1-chome, Shinjuku-ku, Tokyo, Japan
OLYMPUS OPTICAL CO. (EUROPA) GMBH.
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-16, 20097 Hamburg, Germany
OLYMPUS AMERICA INC.
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
OLYMPUS KEYMED
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England
OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373
OLYMPUS BEIJING REPRESENTATIVE OFFICE
Rm. 818, South Tower, Beijing Kerry Centre, No.1 Guanghua Road, Chaoyang District, Beijing, 100020, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
OLYMPUS AUSTRALIA PTY. LTD.
1/104 Ferntree Gully Road, Oakleigh, VIC 3166, Australia

Designed for the toughest cases, this slim scope boasts a user-adjustable flexible insertion tube that's ideal for all colonoscopic applications




innoflex™

Innovative adjustable tube and slim design enable smoother, easier insertion into the colon

Not only does the PCF-160A colonovideoscope feature a slimmer design that reduces patient discomfort, it also incorporates the groundbreaking variable stiffness insertion tube. Thanks to a built-in coil in the insertion tube that can be expanded or contracted as required, the variable stiffness tube is the world's first user-adjustable flexible insertion tube. By allowing you to discretionally modify insertion tube flexibility, the variable stiffness tube reduces the chance of looping and makes insertion beyond the splenic flexure easier. At the same time, it is expected to reduce the need for abdominal pressure on the sigmoid looping and improve access to the cecum.

Range of Stiffness



Adjustment knob for control of tube flexibility

Adjusting the tube's flexibility is quick and easy. All you have to do is turn the knob located just below the control section. Insertion tube flexibility can be adjusted in continuous increments to match the internal conditions and contours of the colon.



Variable stiffness adjustment knob

Anticipated Advantages of Variable Stiffness

- ◎ Improved insertion capability minimizes patient discomfort and reduces the time required for examinations.
- ◎ The cecum can be reached more easily in cases where insertion is difficult.
- ◎ Reduced need for abdominal pressure and position change, reducing staff workload during examinations.

Main Features for Advanced Performance

- The flexibility of the variable stiffness insertion tube can be adjusted to suit the internal conditions and contours of the colon, simply by turning a knob below the control section.
- Slimmer than conventional colonoscopes with a 11.5 mm outer diameter insertion tube, yet still featuring a channel with 3.2 mm inner diameter that's suitable for routine colonoscopy.
- Superior image quality is delivered even in the expanded full height mode when combined with the dedicated CV-160 video system center.
- Extra-wide 140° field of view and extensive 4-way angulation (180° up/down and 160° right/left) enable complete and comprehensive examination of the colon.
- Ergonomically designed grip enhances scope maneuverability while redesigned easy-to-access knobs and user-selectable switches improve operability.
- Fully compatible with the CV-160/140/100.
- Scope ID function stores individual scope information and displays it on the monitor to facilitate endoscopy suite management.



Specifications

Optical System	Field of view	140°
	Direction of view	0° Forward viewing
	Depth of field	3 to 100 mm
Distal End	Outer diameter	11.3 mm
	Outer diameter	11.5 mm
Bending Section	Angulation range	Up 180°, Down 180°, Right 160°, Left 160°
Working Length		L:1680 mm, l:1330 mm
	Total Length	L:2005 mm, l:1655 mm
Instrument Channel	Inner diameter	3.2 mm
	Minimum visible distance	5 mm from the distal end
	Endo-Therapy accessory entrance/exit position in field of view	

