

OLYMPUS[®]


Your Vision, Our Future

Single Use Channel Cleaning Brush

BW-201B/T



Single Use Olympus Cleaning Brushes



As an endoscope manufacturer, Olympus is pleased to offer a brush that has been designed to facilitate maximum cleaning performance while ensuring compatibility with the endoscope channel.

Specifications

	BW-201B	BW-201T
Compatible Channel Diameter	2~3.2mm	2~4.2mm
Working Length	950mm	2200mm

BW-201B/T.A set

Single Use Channel Cleaning Brush

Designed for one-time use, this brush is ideal for cleaning out suction channels in compatible Olympus endoscopes and openings in compatible Endo-Therapy accessories.
Individually packed — 50 packs per box.



BW-201B/T.B set

Channel Cleaning Brush and Channel-Opening Cleaning Brush

This set contains both the Single Use channel cleaning brush and the Single Use channel-opening cleaning brush.
One of each per pack — 50 packs per box.



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS[®]

OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Mondith, 3-1 Nishi-Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914, Japan
OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany
OLYMPUS AMERICA INC.
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
OLYMPUS LATIN AMERICA, INC.
6100 Blue Lagoon Drive, Suite 390 Miami, Florida 33126-2087, U.S.A
KEYMED LTD.
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373
OLYMPUS HONG KONG AND CHINA LIMITED
Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong
OLYMPUS (BEIJING) SALES & SERVICE CO., LTD.
Room 1401, East Tower, Gong Yuan No6 Royal Palace, No6 Gong YuanXijie, JianGuoMenNei, DongCheng District, Beijing, 100005, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
OLYMPUS AUSTRALIA PTY LTD
31 Gilby Road, Mount Waverley, VIC., 3149, Australia

www.olympus.com

Printed in Japan F1111SB-4-0405