

OLYMPUS RETRIEVAL DEVICES

Offering Retrieval Solutions for Any Situation

LEGEND

APPLICATION	PRODUCT
Polyps	FG-45L-1*, 46L-1*, 51D, 52D, 54D, 600U
Food bolus	FG-16L-1*, 51D, 52D, 54D, 600U
Seeds	FG-16L-1*, 51D, 52D
Round, slippery objects	FG-16L-1*, 22Q-1, 23Q-1
Small bones	FG-6L-1*
Flat objects, coins, dentures	FG-4L-1, 9L-1*, 32C-1
Sharp objects (nails, needles, blades, pins)	FG-20P-1*
Pulmonary foreign bodies	FG-26C-1, 51D, 52D, 54D
Grasping tissue in 2-channel EMR, foreign bodies	FG-42L-1*, 50L-1
Soft objects	FG-9L-1*, 14P-1, 26C-1
CBD/PD stones	FG-22Q-1, 23Q-1, 33W, 301Q, 401Q, 402Q, 403Q, V421PR, V431P, V432P, V432PR
CBD/PD Stent removal	FG-22Q-1, 23Q-1, 44NR-1, 402Q, 403Q, V432P, V432PR
CBD/PD Sludge	FG-301Q, 401Q, 402Q, 403Q, V421PR, V431P, V432P, V432PR


* Device also available in colon length.

SINGLE USE

Retrieval Baskets


ERCP

FlowerBasket™ Rotatable




FG-V421PR

FlowerBasket™ Wireguided




FG-V431P

TetraCatch™ Rotatable




FG-V432PR

TetraCatch™ Wireguided




FG-V432P

FlowerBasket™ w/Pin-Vise Handle




FG-401Q

4-Wire Basket




FG-33W

Stiff-Wire Basket



FG-402Q


Soft-Wire Basket



FG-403Q


PULMONARY

3-Wire Spiral Basket



FG-52D


4-Wire Spiral Basket



FG-51D


Grasping Forceps

Polygrab™ Tripod



FG-600U

3-Nail



FG-54D

REUSABLE

Retrieval Baskets

4-Wire Basket w/3-Ring Handle



FG-16L-1

Stiff-Wire Basket w/Pin-Vise Handle



FG-22Q-1

Soft-Wire Basket w/Pin-Vise Handle



FG-23Q-1

FlowerBasket™ w/Pin-Vise Handle



FG-301Q

Grasping Forceps

Rat Tooth Alligator Jaw



FG-44NR-1

Tripod



FG-45L-1

Pentapod



FG-46L-1

Shark Tooth



FG-32C-1

V-Shape



FG-4L-1

Rubber Tip (non-latex)



FG-20P-1

Alligator Jaw



FG-6L-1

Rat Tooth Alligator Jaw



FG-42L-1

Rat Tooth



FG-9L-1

Rat Tooth



FG-14P-1

Rat Tooth



FG-26C-1

Rat Tooth (long and wide)



FG-50L-1

Olympus offers a comprehensive line of devices developed to handle any situation—from extraction of smooth, slippery objects, to removal of sharp, splintered foreign bodies.

Knowledge ➤ Vision ➤ Solutions

To place an order or receive more information, please contact customer service at 800-848-9024 or visit olympusamerica.com/endothorapy