

# OLYMPUS®



EVIS EXERA GASTROINTESTINAL VIDEOSCOPE  
**OLYMPUS GIF TYPE Q160**


# LARGE DISPLAY WITH SUPERIOR IMAGE QUALITY

Large images with superior sharpness and clarity, expanded field of view, and great depth of field make this high-quality scope the right choice for a wide variety of gastroscopic examination and treatment procedures





## Main Features for Superb Performance

- Large display size delivers superior image quality with sharper, crisper image and true-to-life color from edge to edge.
- Extra-wide 140° field of view enables accurate observation of a wider area.
- Slim 9.5 mm outer diameter insertion tube with 2.8 mm channel diameter.
- 4-way angulation (210° up, 90° down and 100° right/left) enables complete and comprehensive examination of the upper digestive tract.
- Ergonomically designed grip enhances scope maneuverability while redesigned easy-to-access knobs and user-selectable switches improve operability.
- Fully compatible with the CV-160/140.
- Scope ID function stores individual scope information and displays it on the monitor to facilitate endoscopy suite management.





### Specifications

Optical System	Field of view	140 <sub>i</sub>
	Direction of view	0 <sub>j</sub> Forward viewing
	Depth of field	3 to 100 mm
Distal End	Outer diameter	9.8 mm
	Insertion Tube	Outer diameter
Bending Section	Angulation range	Up 210 <sub>i</sub> , Down 90 <sub>i</sub> , Right 100 <sub>i</sub> , Left 100 <sub>i</sub>
Working Length		1030 mm
Total Length		1345 mm
Instrument Channel	Inner diameter	2.8 mm
	Minimum visible distance	3 mm from the distal end
	Endo-Therapy accessory entrance/exit position in field of view 	



### Olympus business areas



**OLYMPUS**  
FOCUS ON LIFE

<http://www.olympus.com>

#### OLYMPUS OPTICAL CO., LTD.

San-Ei Building, 22-2, Nishi Shinjuku 1-chome, Shinjuku-ku, Tokyo, Japan

OLYMPUS OPTICAL CO. (EUROPA) GMBH,  
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-16, 20097 Hamburg, Germany

OLYMPUS AMERICA INC.,  
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.

OLYMPUS KEYMED  
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.  
491B River Valley Road #12-01/04 Valley Point Office Tower, Singapore 248373

OLYMPUS BEIJING REPRESENTATIVE OFFICE  
Rm. 818, South Tower, Beijing Kerry Centre, No.1 Guanghua Road, Chaoyang District, Beijing, 100020, China

OLYMPUS MOSCOW LIMITED LIABILITY COMPANY  
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia

OLYMPUS AUSTRALIA PTY. LTD.  
1/104 Ferntree Gully Road, Oakleigh, VIC 3166, Australia