

OLYMPUS[®]

Your Vision, Our Future



OES Bronchofiberscope
OLYMPUS BF TYPE 1T60

Maximum therapeutic and suction capabilities
with an extra-large channel




OES Bronchofiberscope

OLYMPUS BF TYPE 1T60




Extra-wide channel accommodates 2.8 mm diameter forceps, greatly enhancing therapeutic capability



- Extra-wide instrument channel measuring 3.0 mm across enables use of 2.8 mm diameter forceps and increases suction power, while retaining the slim 5.9 mm distal end diameter of its predecessor.
- Wide 120° field of view enables accurate observation of a wider area.
- 180° up/130° down angulation for easier maneuvering in the bronchi.
- 600 mm working length ideal for use in the tracheobronchial tree.
- High-frequency electrosurgical devices can be used with this endoscope.
- Compatible with both OES and EVIS video systems.

■ Specifications

Optical System	Field of view	120°
	Direction of view	0° (Forward viewing)
Insertion Tube	Depth of insertion	~ 600mm
	Distal end outer diameter	~ 5.9mm
	Insertion tube outer diameter	6.0mm
Instrument Channel	Working length	600mm
	Channel inner diameter	3.0mm
Bending Section	Minimum visible distance	5mm from distal end
	Angulation range	Up180°, Down130°



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS®

OLYMPUS CORPORATION
Shinjuku Monolith, 3-1 Nishi-Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914, Japan
OLYMPUS WINTER & IBE GMBH
Kuehnstrasse 61, 22045 Hamburg, Germany
OLYMPUS EUROPA GMBH
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany
OLYMPUS AMERICA INC.
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
OLYMPUS LATIN AMERICA, INC.
6100 Blue Lagoon Drive, Suite 390 Miami, Florida 33128-2087, U.S.A
KEYMED LTD.
Jeymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.
831B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248072
OLYMPUS HONG KONG AND CHINA LIMITED
Room 1820-1821, Ocean Centre, 3 Canton Road, Tsimshatsui, Kowloon, Hong Kong
OLYMPUS BEIJING REPRESENTATIVE OFFICE
Room 1401, East Tower, Gong Yuan Hotel Royal Palace, Noli Gon YuanXike, JianGuoMenWai, DongCheng District, Beijing, 100006, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kuznetzkaya 13, bld. 1, B-2, Russia
OLYMPUS AUSTRALIA PTY LTD.
31 Gilby Road, Mount Waverley, VIC, 3149, Australia

www.olympus.com

Printed in Japan R395SB-3-0604