

TC-EU COMPACT ENDOSCOPIC ULTRASOUND WORKSTATION

A DEDICATED COMPACT ENDOSCOPIC ULTRASOUND WORKSTATION

The Olympus Endoscopy Support range of products includes a complete family of workstations, developed to be practical and space-efficient wherever they are used. Ergonomically designed, with many innovative features, they are more than just trolleys, they are unique, purpose-designed mobile workstations.

The Olympus TC-EU is a compact workstation designed to support the EU-M30S Endoscopic Ultrasound processor.

FEATURES INCLUDE:

- **Compact** - space-efficient for use in what is so often a cramped working environment.
- **Aesthetically pleasing design** - to suit the modern clinical environment.
- **Base tray** - designed to accommodate ancillary items for the EU-M30S.
- **Mobile** - with four 100mm castors, two diagonally opposed with brakes. Easily brought into or moved out of the procedure area.
- **Colour co-ordinated bumper protection on all four corners** - to prevent accidental damage to walls, doors and paintwork.
- **Keyboard shelf** - easy keyboard access for operation of EU-M30S processor.
- **Top shelf** - to accommodate a 14" monitor.
- **Fixed second shelf** - to accommodate the EU-M30S processor.
- **Adjustable third shelf** - to accommodate user specified equipment, such as video printer or video recorder.
- **Integral handle** - ideally positioned for easy manoeuvrability.
- **Back Panel** - is available to provide a fully enclosed workstation.
- **Plug Blocks** - a multiway plug block kit is available reducing the number of trailing cables

CONSTRUCTION:

Olympus workstations, designed and manufactured by KeyMed, have an enviable reputation for performance and durability the TC-EU is no exception.

Built to an exhaustive quality standard, the all metal frame and shelf construction are electrostatically powder painted and will be highly durable within the operating environment.



SPECIFICATION:

Dimensions: 480(W) x 1150(H) x 570(D) mm

Unladen Weight: 29 kg

Maximum carrying capacity: 100 kg

Olympus reserves the right to alter the specification without notice.

OLYMPUS®

OLYMPUS OPTICAL CO., LTD.

San-Ei Building, 22-2, Nishi Shinjuku 1-chome, Shinjuku-ku, Tokyo, Japan

OLYMPUS KEYMED

KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom

OLYMPUS OPTICAL CO., (EUROPA) GMBH.

(Premises/Goods delivery) Wendenstraße 14-18, 20097 Hamburg, Germany
(Letters) Postfach 10 49 08, D-20034 Hamburg, Germany

OLYMPUS AMERICA INC.

Two Corporate Center Drive, Melville, New York 11747-3157, USA

OLYMPUS SINGAPORE PTE LTD.

491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore



ISO 9001 Registered
Certificate No. FM 20993

OLYMPUS

TC-EU COMPACT ENDOSCOPIC ULTRASOUND WORKSTATION



**A further product in a design
co-ordinated endoscopy support range**

